



Vertical Glass Distillation Quartz Heater

Vertical Glass Distillation Quartz Heater



ACMAS Technologies Inc.

CE An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603
Email: acmastechnologies@gmail.com | obromax.delhi@gmail.com
www.acmastechnologies.com | www.acmas.in

ACMAS Technologies Inc.

CE An ISO 9001:2000 Company

These units consist of main boiler. A built in heater enclosed in Quartz glass tube. Housing auto Water level adjuster, Bottom outlet for cleaning waste discharge which is fitted with Teflon screw Cap valve. A coil condenser fitted on top for condensing vapours. The whole unit is mounted on epoxy coated metal box stand & covered with acrylic sheet cover for protection. Power supply unit is in built no extra attachment is required. Housing sturdy & reliable power System require 1 K.W. Power per heater.

Vertical Glass Distillation Quartz Heater

Vertical Glass Distillation Quartz Heater

Features:

- Temperature resistant
- Saving of electrical energy in every start
- Convenient bottom drain for easy drainage as accumulated solids impair the efficiency, cleaning is also possible

- Distilled water in cold condition
- Automatic start, cut-off & re-start in synchrony with water pressure/opening by electronic circuit.

