



**Rectangular Hot Plate
With Cast iron**

Rectangular Hot Plate With Cast iron

www.acmastechnologies.com | www.acmas.in



ACMAS Technologies Inc.

CE | An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603
Email: acmastechnologies@gmail.com | obromax.delhi@gmail.com
www.acmastechnologies.com | www.acmas.in

ACMAS Technologies Inc.

CE | An ISO 9001:2000 Company



Rectangular Hot Plate With Cast iron

www.acmastechnologies.com | www.acmas.in



Top plate made from S.S. Heavy Sturdy Cast Iron Casting and coated with epoxy base resistive coating that does not chip off as ordinary paint even at higher temperatures.

Maximum temperatures with thermostat 275°C and 350°C with energy regulator.

Rectangular Hot Plate With Cast iron

www.acmastechnologies.com | www.acmas.in

DESCRIPTION

- Rectangular shape with cast iron or stainless steel top.
- This is a heavy duty heating plate with a removable cast iron plate / Stainless steel Plate at the top.
- The housing is made of mild steel sheet and finished in powder coat paint

TECHNICAL DATA

- A heat control switch for low, medium and high heat or energy regulator and two nos. indicating lamps are provided on the front panel

- The hot plate has a maximum power consumption of 3.0 K.W.
- The unit is designed for operating on 230 V. Single phase 50 C/S. mains.

APPLICATION

- Medical Lab, Testing Laboratory, Pharmaceutical industries etc.

