



VDRL Rotator



VDRL Rotator



OBROMAX

ACMAS Technologies Inc.

CE An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603
Email: acmastechnologies@gmail.com | obromax.delhi@gmail.com
www.acmastechnologies.com | www.acmas.in



ACMAS Technologies Inc.

CE An ISO 9001:2000 Company



OBROMAX

VDRL Rotator

Our manufactured and export VDRL Rotator is suitable for VDRL tests, blood grouping tests and for mixing solutions in small flasks, breaker or bottles. Made out of finest raw material available, it provides noiseless and smooth movement. Tray type platform size 300x300 mm, with great motor durability and speed control switch. Due to our commitment in meeting the demands of our clients, we have emerged as the most trusted exporters and suppliers of this product. Our priority is to satisfy the clients' needs.

- Platform : 300x300mm
- Speed : Control Switch
- RPM range – upto 250 RPM
- Electronic control
- 0-60 min timer.
- Workable on 220 V AC 50 Hz Single phase.



VDRL Rotator