

Water still wall mounting type

Wall Mounted

Water still wall mounting type

Wall Mounted



ACMAS Technologies Inc.

CE An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603
Email: acmastechnologies@gmail.com | obromax.delhi@gmail.com
www.acmastechnologies.com | www.acmas.in

ACMAS Technologies Inc.

CE An ISO 9001:2000 Company

Water still wall mounting type

Water still wall mounting type

This unit is recommended for making pyrogen free water of a very high degree of purity.

Operation

The water still designed for wall mounting and have an identical method of operation. Water feeds at the base through tap of the condenser jacket. Water circulates around the condenser tube and after that fills the boiling chamber. The height of the water is controlled by the overflow level. After this process switch on the heaters. For efficient operation, a quantity of cooling water is required equal to approximately eight times of the distillate output.



Construction

- Water still is made of stainless steel sheet.
- The portable condenser system is provided for easy cleanliness fo inside portion for better efficiency & avoids choking.
- 7 Electrically heated water still are preferable to gas heated waterstill for airconditioned laboratories.

Electric supply

To be connected only a socket of 230 V AC 50 Hz supply.

Capacity	Load
2 ltrs per hr	1.5 KW
4 ltrs per hr	3.0 KW
6 ltrs per hr	4.0 KW
8 ltrs per hr	6.0 KW
10 ltrs per hr	8.0 KW