



# Laboratory Stirrer



**ACMAS Technologies Inc.**

CE | An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA  
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603  
Email: [acmastechnologies@gmail.com](mailto:acmastechnologies@gmail.com) | [obromax.delhi@gmail.com](mailto:obromax.delhi@gmail.com)  
[www.acmastechnologies.com](http://www.acmastechnologies.com) | [www.acmas.in](http://www.acmas.in)

**ACMAS Technologies Inc.**

CE | An ISO 9001:2000 Company



# Laboratory Stirrer



## Laboratory Stirrer

Our entire Laboratory Stirrer Range is considered as congregate requirements of laboratories challenging wide elasticity of stirring operations. Our Electric Lab Stirrer has a chuck to enable quick interchangeability of shafts. There is different DIA 1/20hp up to 6mm. All our models are with 1/8 hp AC/DC motor and chuck for speedy shaft interchangeability up to 8 mm Dia. One of the reputed Lab Stirrer Manufacturers and Suppliers we have made it a point to ensure our customers of our optimum quality and high performance of our products.

### Properties

- Has a chuck to enable quick interchangeability of shafts
- There is different DIA 1/20hp up to 6mm
- Our models come with 1/8 hp AC/DC motor and chuck for speedy shaft interchangeability up to 8 mm diameter

