



## Reciprocating Khan Shaking Machine



**OBROMAX**

**ACMAS Technologies Inc.**



CE An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA  
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603  
Email: [acmastechnologies@gmail.com](mailto:acmastechnologies@gmail.com) | [obromax.delhi@gmail.com](mailto:obromax.delhi@gmail.com)  
[www.acmastechnologies.com](http://www.acmastechnologies.com) | [www.acmas.in](http://www.acmas.in)

**ACMAS Technologies Inc.**

CE An ISO 9001:2000 Company

**OBROMAX**



## Reciprocating Khan Shaking Machine

## Reciprocating Khan Shaking Machine

### Reciprocating Khan Shaking Machine

Kahn racks, bottles or flasks of different sizes can be shaken by adjusting the compartments provided on the shaking rack. Minor rate of motion can be adjusted mechanically by belt tension from 150 to 200 r.p.m. approximately. Fitted with electrically operated motor to work on 220/230 Volts A.C. With Perfect / Crompton Motor

Plant form size	Body	Speed Motor Capacity
(a) 20" x 13" x 4"	Closed	Fixed 1/4 H.P
(b) 31" x 17" x 4"	Closed	Variable 1/4 H.P
(c) 36" x 18" x 4"	Closed	Variable 1/2 H.P

