



ACMAS Technologies Inc.

CE | An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603
Email: acmastechnologies@gmail.com | obromax.delhi@gmail.com
www.acmastechnologies.com | www.acmas.in

Wrist action shaking machine



ACMAS Technologies Inc.

CE | An ISO 9001:2000 Company





www.acmastechnologies.com | www.acmas.in

platform sizes	holding capacity	h. P. Motor
45 x 45 cms.	25 flasks of 100 ml.	0.25
45 x 45 cms.	36 flasks of 50 ml.	0.25
45 x 45 cms.	16 flasks of 250 ml.	0.25
55 x 55 cms.	25 flasks of 250 ml.	0.25
65 x 65 cms.	36 flasks of 250 ml.	0.50
100 x 100 cms.	81 flasks of 250 ml.	1.00
100 x 100 cms.	81 flasks of 500 ml.	1.00

Wrist action shaking machine

Continuous shaking of flasks from 100 to 500 ml capacity. Bench type table top model with platform with rubber discs/steel clamps mounted on ball-bearing crank-shafts to give an orbital motion to flasks within a diameter of 2". Motor is coupled with v-belt with a pitch pulley and mounted on heavy angle iron frame. Machine supplied with fixed or variable speed controlled by mechanical speed controller. To work on 220 v ac 50 hz single phase.



www.acmastechnologies.com | www.acmas.in