



Mortuary Chamber



ACMAS Technologies Inc.

CE An ISO 9001:2000 Company





Introduction

Obromax Mortuary Chambers are designed for storing cadaverous under cool or low temperature conditions to prevent decomposition while preserve in Hospitals, Morgue, Railway & Airport. These chambers offer vast space savings advantages over the conventional type. These chambers also maintain greater hygiene. Our Mortuary Chambers are highly appreciated among our customers & served with best pre & post technical services.

Salient features:

- Customized storage (1/2/3/4/5-as per customer's requirement)
- CFC free & Eco-friendly refrigerants
- Equipped with coaxial blowers, quality compressors & evaporators
- Double walled chamber
- Inner body made up of high grade stainless steel
- Outer body made up of epoxy powder coated mild steel or can be customized
- Simplified cadaver handling
- Prevent decomposition
- Insulated door
- Temperature range 2 to 8 deg C/ -15 to 10 deg C or as per customer's requirement.
- Operating voltage: 220 V, 50 Hz, single phase
- GMP & GLP compliant
- CE certified
- Trademark
- Factory calibrated certificates

Optional Accessories:

- Timer (999 min) • Real Time programmer • Data Acquisition Program • Serial Data Port
- Inspection Window in Door • Printer report program • Access Port

Technical Matrix

Construction	Inner	SS-304
	Outer	Powder Coated MS
	Doors	Inner wall is of SS 304 & outer is of powder coated MS
Temperature	Castors	Lockable
	Controller	PID or PLC (as required)
	Sensor	PT 100 or PT 1000 (as required)
	Range	2 to 8 deg C or -15 to 10 deg C (Customizable)
	Dsisplay	LCD/LED
Illumination System		Door operated illumination system comprises of Fluorescent lights
Adj. Alarm limits		Optional
Condenser	AirCirculation	Efficient condenser with automatic evaporating system (condensate) Forced
	Refrigerant	CFC free or R134
Water reservoir		18 ltrs (Customizable)
Light Intensity		As per requirement
Capacity		1/2/3/4/5/6 (Customizable)
Internal dimensions		Model Specific
Maximum Permitted Load per shelves		155 kg
Power Supply		220 V & 50/60 Hz

ACMAS Technologies Inc.

 An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA

Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603 | Email: acmastechologies@gmail.com | obromax.delhi@gmail.com
www.acmastechologies.com | www.acmas.in